

# GEMINO




**GEMINO 30** walker • **GEMINO 30M** walker

**Rollators providing excellent support and stability  
– ideal for walking training.**

# GEMINO 30 walker

# GEMINO 30 M walker

Gemino walker is almost identical to the standard Gemino 30 model, but provides extra support and stability.



Drip holder



Crutch holder



Light



Table



Umbrella



Slow-down brake



Bag for oxygen bottle



Bag with lid and shoulder strap

## Comfortable, correct position for walking

The unique and very stable design ensures that you can walk securely in all settings.

The forearm supports can be easily adjusted for height, position and angle to suit the individual user.

The brake and pushing handles are integrated into the forearm supports and are designed to be comfortable and easy to use.

## Safe and easy to use

High-quality integrated reflectors and a reflective band on the bag make Gemino clearly visible to road users.

The unique brake block design ensures that the brakes stop the wheels very effectively.

The design of the tipping bar makes it very easy to negotiate obstacles such as kerbs and door sills. A unique side deflector makes it easier to navigate doorways, table legs and other furniture, without the danger of the back wheels getting caught.

## Gemino – design

The items with which we surround ourselves help define how we are perceived. Gemino's elegant design appeals to both young and old alike.

## Wide range of accessories

Gemino has a wide range of accessories. Together with the adaptation options in the rollators themselves, these accessories make it possible to tailor Gemino to match user requirements exactly.

Specifications	Gemino 30 walker	Gemino 30 M walker
Total height	98–115 cm	84–100 cm
Total length	65 cm	65 cm
Total width	61 cm	61 cm
Seat height	62 cm	55 cm
Total weight	10.9 kg	10.7 kg
Max. user weight	150 kg	130 kg
Max. weight in bag	5 kg	5 kg
Recommended user height	150–200 cm	135–170 cm



handicare

[www.handicare.com](http://www.handicare.com)

For complete product information, visit your local dealer or contact customer service:

T +31 (0)492 593 888  
F +31 (0)492 537 931  
[info@handicare.com](mailto:info@handicare.com)