



**PRO** **ACTIV<sup>®</sup>**

**LIFT** **activ**  
mechanical

INCREDIBLY FLEXIBLE AT WORK & IN EVERYDAY LIFE

# LIFT activ

The **LIFT activ** is equipped with a manual seat height adjustment option, which is easily adjusted using an ergonomic handle. The continuous seat height adjustment is achieved with gas pressure springs. With only a little muscle power, the seat is lifted upwards by up to 30 cm. In this position, the wheelchair user is eye to eye with a person standing about 1.75 m tall.

This adjustable seat height makes it highly flexible in use in a wide range of professions, and lets you use conventional furniture and kitchen facilities. With the LIFT activ, you can work flexibly at different workstations, machines and devices.



## LIGHT & ACTIVE DESPITE THE LIFTING MECHANISM

The frame and lifting mechanism of the LIFT activ impresses with a weight starting at 14.5 kg and with an attractive and dynamic design.



## SIMPLE TRANSPORT

Transport in a car boot is problem-free - simply take the wheelchair apart into its individual modules, with no need for tools.



## WIDE RANGE OF OPTIONAL EXTRAS & ACCESSORIES

Each LIFT activ is a customised model, with frame dimensions and equipment that are always tailored to the wheelchair user's needs. You can choose from a wide range of equipment for the LIFT activ, including footrests, brakes, push handles, back and seat systems, as well as drive and caster wheels.



## THE LIFT ACTIV GETS YOU TO THE RIGHT LEVEL - AT WORK, IN EVERYDAY LIFE & LEISURE

Thanks to its low weight (in comparison to the LIFT solid), the **LIFT activ** is ideal for transport and for using the wheelchair at work, in everyday life and leisure activities. Reaching the files on the top shelf at the office in the morning, putting the dishes in the wall cupboard above the worktop in the kitchen at home at lunchtime, and having a conversation with family and friends at eye level at the bar in the evening – that's no problem with the LIFT activ.



## FLEXIBLE OPTIONS FOR USE - USE AS A STANDING-UP AID

For people whose conditions leave them some movement in their legs, the LIFT activ can also be used as a standing-up aid.

For this use, it makes sense to select the “Divided footrest with swivel frame stem” or another footrest that can be folded or swivelled away.

### OPTION: REMOVABLE FRAME STEM, FOLDABLE WITH DIVIDED FOOTREST

The frame stem can also be swivelled out and removed, so that the frame ends at the front end of the cover and leaves a lot of space for transfer.



## TECHNICAL SPECIFICATIONS

Lower leg guidance:	G- or V-shape frame stems available
Front frame angle:	100°   105°
Lift function/seat height adjustment:	Manual seat height adjustment, over 30 cm infinitely adjustable
Seat width:	35 - 46 cm (in 1 cm steps)
Seat depth:	35 - 48 cm (in 1 cm steps)
Back height:	20 - 48 cm (in 1 cm steps)
Lower leg length:	30 - 50 cm (in 1 cm steps)
Rear seat height:	41 - 43 cm (in 1 cm steps)
Front seat height:	45 - 51 cm (in 1 cm steps)
Tipping point:	6 - 10 cm (in 1 cm steps)
Total product weight:	Starting at 14.5 kg
Max. standard load weight:	Up to 100 kg
Possible drive wheel camber:	1°, 4°, 2.5°