

JUVO care™

Hygiene

Mobile shower commodes & accessories







Tilt and recline

5000 series



Tilt-in-space Reclining backrest

Seat height adjustable, configurable tilt-in-space



rest version



Adjustable seat height
6 settings at 26 mm



Fully adjustable
seat tilt







Removable backrest for
headrest option

Specifications

SKU	Standard		Bariatric	
	JH5041ATAR	JH5046ATAR	JH5051ATAR	JH5061ATAR
Seat width	410	460	510	610
Width between arms	440	490	540	610
Seat height anterior 6 settings at 26mm	510–660			
Seat height posterior 6 settings at 26mm	480–630			
Seat to heel	320×575			
Seat tilt angle	-5° to +38°			
Width overall	560		610	730
Height overall	1465–1615			
Depth overall — tilted	940–1600			
Armrest height from seat, low-high	200–240			
Backrest height and angle adjustment	800mm from seat and 140 degree recline from seat			
Aperture, wall to rear	300			
Toilet bowl clearance height	450–600			
Toilet bowl clearance width	400	410	470	500
Toilet bowl clearance between pan carrier	315			
Max user weight, kg	200		300	
Chair weight excluding seat, kg	24	25	26	28

Commode dimension reference guide



Materials & warranty

Part	Warranty period, years	Material	Grade	Finish
Frame	7	Stainless steel	304	Mirror
Castors	2	Stainless steel	304	Mirror
Fixings, bolts, nuts	2	Stainless steel	304/316	Mirror
Spring & button	2	Stainless steel	304/316	Mirror
Seat	2	Polyurethane	Anti-Bacterial	Smooth
Backrest straps	2	Polyurethane	Anti-Bacterial	Smooth
Hand grips	2	Polyurethane	Anti-Bacterial	Non-slip
Rear Wheels, self-propelled	2	Moulded plastic	Anti-Bacterial	Smooth
Armrests	2	Polyurethane	Anti-Bacterial	Smooth



Always innovating...

Juvo Solutions PTY LTD

ABN 93 151 426 416

+61 8 6454 9960

enquiries@juvo.com.au

www.juvo.com.au



Australian
Standard

AS/NZS 3973:2009

TGA

Juvo Care
Hygiene products
are entered in
the ARTG, 191619

Authorised Distributor